



designed to work perfectly

German Technology made in Brazil

/// presentado por Natalia Casadiego



/// IKA Brasil

IKA

designed to work perfectly

Campinas- São Paulo- Brazil



IKA

designed to work perfectly

Lab dancer

/// Test Tube Shaker



TOUCH FUNCTION!

Attractively designed, test tube shaker.

Designed for mixing small test samples with a touch function.

- Small, compact and reliable
- For small containers up to 30 mm in diameter, e.g. test tubes, centrifuge tubes, Eppendorf beakers
- **Fixed speed (2800 rpm)**
- Excellent mixing action

Vortex 3

///Vortex shaker suitable for short-time operation



Different applications thanks to 3 interchangeable attachments and 7 inserts: Eppendorf tubes, microtiter plates, Erlenmeyer flasks 250 ml etc...



- Shaking movement: orbital
- Orbital diameter: 4 mm
- Infinitely adjustable speed range: 500 - 2.500 rpm






IKA

designed to work perfectly

Attention!

Falcon 15ml – 50 ml



Designation	Description	Id.- No.	Touch Mode	Continuous Mode	Range of speed
 V 3.1 Standard attachment	For test tubes and small vessels upto diameter 30 mm	33 412 00	X	X	
 V 3.2 One- hand attachment	88 mm-plate attachment with rubber pad	33 423 00	X	X	
 V 3.3 Universal attachment	150 mm-plate attachment with rubber pad	33 424 00	-	X	



Petri dish

MS 3 Digital

/// Compact, universal small shaker suitable for shaking tasks with all small vessels and microtiter plates



- Speed Range: 0- 3000 rpm (mode A, B)
- Shaker Diameter: 4.5 mm
- Timer with countdown function
- Continuous or touch operation (with standard attachment)

Scope of Delivery:



MS 3.1 Standard attachment



MS 3.3 Universal attachment



MS 3.4 Microtiter attachment



MS 1.21 One-hand Insert



MS 1.32 Test tube Insert




designed to work perfectly

MS 3 Digital

Customer want to ensure, whether we could use MS 3.5 to shake, the PCR plates would not fly out



MS 3.5



KS 501 Digital

/// Low profile laboratory orbital shaker with large mounting surface and load capacity of up to 15 kg.



AS 501.1




- Infinitely variable speed control from 0-300 rpm
- Digital display
- Orbital diameter 30 mm
- Ideal for vessels with a volume of more than 250 ml, e.g. round flasks, Erlenmeyer flasks, culture flasks and bottles
- Ideal for biological and microbiological growth tests and for production of stock solutions
- Timer function
- **Integrated USB and RS 232 interface!!! → New!**

IKA

designed to work perfectly



10 ml de Penicilin
KS 130 Digital



Test Tube
KS 130 Digital



Test Tube 30 ml
KS 130 Digital



20 ml Bottle
KS 130 Digital



8 bottles – 250 ml/ Filled with 200 ml
KS 260

Magnetic Stirrers

Resources



Lab disc / Color
squid / Topolino



Pequenos volumes
Sem aquecimento
Sem display

RH basic/digital



ext. Sensor: **ETS-D**
Interface: **no**
Timer / Counter: **no**
Display: **LED**

C-MAG basic



ext. Sensor: **ETS-D**
Interface: **no**
Display: **LED**

RCT basic



ext. Sensor **PT1000** / **ETS-D**
Interface: **não**
Timer / Counter: **no** / **no**
Display: **LED**

RET basic



ext. Sensor **PT1000** / **ETS-D**
Interface: **não**
Timer / Counter: **no** / **no**
Display: **LED**

C-MAG digital



ext. Sensor: **PT1000**
Interface: **no**
Display: **LCD**

RET control-visc



ext. Sensor **PT100**
Interface: **USB** / **RS232**
Timer / Counter: **Si** / **S**
Display: **TFT**

C-MAG HS 7 Control



ext. Sensor **PT1000**
Interface: **USB** / **RS232**
Timer / Counter: **Si** / **Si**
Display: **LCD**

RCT Digital



ext. Sensor **PT1000**
Interface: **USB** / **RS232**
Timer / Counter: **Si** / **Si**
Display: **LCD**

Price



Color Squid

///Small magnetic stirrer with glass surface



- Digital speed display (LED)→ the client can see the current rpm
- Electronically controlled motor for more capacity
- Higher speed range from 0 - 2.500 rpm
- Max. stirring quantity 1 l
- Stirring quantity max. per stirring position (H₂O): 1 L
- Speed range: 0 - 2500 rpm

C-MAG HS Series

/// Magnetic stirrer with heating and ceramic heating plate

Technical Data	HS 4	HS 7 	HS 10 
Temperature	50 - 500 °C	50 - 500 °C	50 - 500 °C
Set-up plate material	Vidro-cerâmica	Vidro-cerâmica	Vidro-cerâmica
Set-up plate dimensions			
Stirring quantity max. per stirring position (H ₂ O)	5 Lts	10 Lts	15 Lts
Speed range	100 – 1500 rpm	100 – 1500 rpm	100 – 1500 rpm
Connection for ext. temperature sensor	no	ETS-D5	ETS-D5





Item No.	Qty.	Description
	1	IKA C-MAG HS 4 / HS7 / HS10
3378000		ETS D5
2437700	1	H 44
1545100	1	H 16V
3547700	1	H 38



ETS-D5



H 38 Holding rod



H 44 Boss head clamp



H 16 V Support rod



designed to work perfectly

ETS-D5

///Electronic contact thermometer

- Incl. stainless steel sensor H 62
- Recommended with :RCT Basic Safety Control, RET Basic Safety Control and RET Basic Safety Control- DIN 12878, class 2
- Temperature Measurement Range: - 50°C -450 °C
- Accuracy of temperature measurement ± 0.2 + tolerance PT1000 (DIN IEC 751 Class A) K

- Advantages:

there is no relevant technical difference between the two devices; the HS7 digital is actually an HS7 basic with the ETS-D5 built into the device, same hardware, same software. The devices should therefore not differ in terms of control behavior for similar applications.



IKA



designed to work perfectly

Coming soon!

CMAG HS Digital

\\Hotplate made of ceramic glass



- HS 4 Digital
- HS 7 Digital 
- HS 10 Digital 
- Heating Temperature range: 50-500°C
- Directly connection for a PT 1000 temperature sensor enable a precise temperature control (included in delivery)
- Control accuracy in medium +/- 0,5 K (in combination with PT 1000)



designed to work perfectly

Comparing C-MAG HS with CMAG HS 7 Digital

[Copy of HS7 Basic VS HS7 Digital Performance Exp test results.xlsx](#)

- Advantages:

there is no relevant technical difference between the two devices; the HS7 digital is actually an HS7 basic with the ETS-D5 built into the device, same hardware, same software. The devices should therefore not differ in terms of control behavior for similar applications.

1. promote **CMAG HS 7 Basic (only main unit)** according to the budget→ ETS-D5 will upgrade your device with direct and accurate control of sample temperature



IKA

designed to work perfectly

CMAG HP Series

\\Hotplate made of ceramic glass



- HP 4
- HP 7 
- HP 10 
- Heating Temperature range: 50-500°C



C-MAG HP 7 + ETS-D5+
Eurostar 400 Control
preparing Polymer



designed to work perfectly

Application

The cliente want to do the red wine distillation test with the HP7



Heating time	Initial Volume	Final Volume
45 minutes	150 ml	62 ml

Sample preparation:

100 ml → Ethanol 12%

50 ml → Agua

RCT Basic

/// The bestseller in laboratories: Strong motor for a higher range of speeds, additional temperature control mode for faster heating of medium

Speed range	50 - 1500 rpm
Heating temperature range	Room temp. + device self heating - 310 °C
Set-up plate material	Aluminium alloy
Set-up plate dimensions	Ø 135 mm
Stirring quantity max. per stirring position (H2O)	20 Litros
Heat output	600 W



**Included on delivery PT 1000.60+
H 100**

Tips: **RET Basic** Stainless steel set-up



PT 1000.60 Temperature sensor, stainless steel





designed to work perfectly


5-10-15 Position Magnetic Stirrers

///Digital magnetic hotplate stirrer, designed for synchronous heating and stirring

Versions RT 5 / RT 10 / RT 15 → Volume/Position: 0,4 Lts



The magnetic coil technology provides noiseless and consistent stirring on all positions. The surface temperature is infinitely adjustable up to 120 °C, producing a maximum medium temperature of 70 °C (depending on the used vessel). Speed remains constant, even when load changes



Overhead Stirrers



Mixing		Disperse phase		
		Solid	Liquid	Gas
Continuous phase	Solid	Solid mixing	Extrusion	Fluidisation
	Liquid	Suspension	Emulsification	Aeration
	Gas	Swirling	Atomization	Gas mixture (homogeneous)

I K A O v e r h e a d S t i r r e r s



Disperse Systems

- The single phases are not or hardly dissolved in each other and separated by a phase boundary

Continuous phase

Disperse phase





designed to work perfectly

RW 20 digital

/// The bestseller in the laboratory










Rpm 1	60 - 500 rpm
Rpm 2	240 – 2000 rpm
Stirring quantity max. per stirring position (H2O)	20 litros
Speed display	Digital
Viscosity max.	10,000 cps
Potencia do motor	70/35 W
Torque max. at stirring shaft	150 ncm



<u>Part Number</u>	<u>Package</u>
3593000	RW20 Digital
3160100	R1826 Stand
2657700	R182
3008600	RH3



Elements Stirrers

Stirrer typ	Design	flow			task					speed		
		axial	radial	tangential	homogenization	Suspension	Emulsification	aeration	heat exchange	low	medium	high
Propeller stirrer		x			x	x			x		x	x
Centrifugal stirrer		x			x	x			x		x	x
Dissolver stirrer			x				x	x			x	x
Paddle stirrer				x	x	x		x		x	x	
Anchor stirrer				x					x	x		
spiral stirrer		x			x				x	x	x	
Blade stirrer			x		x		x	x	x		x	x
Turbine stirrer		x				x	x			x	x	x
Moebius stirrer		x				x	x			x	x	





designed to work perfectly

Eurostar 60 Digital/ Control

/// Universally used strong laboratory stirrer for quantities up to 40 l (H₂O)

Stirring quantity max. (H₂O): 40 l

Max. Viscosity: 50,000 mPas

Speed range: 0/30 – 2000 rpm

Max. torque at stirring shaft: 60 Ncm

Digital vs Control:

- LED vs TFT
- Intermittent operation only control version
- Temperature Sensor: PT 1000 only control, versio..
- Torque Trend, Timer, Temperature Measurement: only control



IKA

designed to work perfectly

Application



70L electronic slurry,
using Eurostar 200
control and R1342
propeller, 1200RPM
continuous operation
for 1 week.

Grinding & Dispersing





designed to work perfectly

Analytical Mill

/// Batch mill

Impact grinding of hard, brittle or non-elastic grinding materials with high-grade stainless steel beater. This beater can be used for a Mohs hardness up to 6 .

Cutting grinding for pulverizing soft, fibrous materials with a cutting blade

Item No.	Qty	Description
2900000	1	IKA. Including: (1) A11.1 stainless steel (1) A11.5 chamber 80 ml
		Accessories :
2904600	1	A11.1 spare beater(hasta 5 Mohs)
2905200	1	A11.2 cutting blade (fibrous grinding materials)
2983000	1	A11.3 beater(9 Mohs)
2904100	1	A11.4 chamber 250ml (Acero Inox. 316 Makrolon)
2983100	1	A11.5 chamber 100ml (Acero Inox. 316 Makrolon)



Speed max: 25000 rpm
Speed max: 80 ml
Power-ON time: 1 min
Power-OFF time: 10 min

Tube Mill control

/// Patented!

Larger cover for larger grinding chamber
From 40 ml to 100 ml

Less cleaning costs
Batch mill with disposable grinding chambers
Eliminate possibility of cross-contamination
Eliminate the effort of cleaning

Continuous monitoring of grinding tests
Thanks to transparent chamber and cover

USB interface
Control and document all parameters





designed to work perfectly

TUBE MILL 100 control

Accessories

Included on delivery scope

Tube			Feature	Benefit
Disposable tube	● MT 40.10 MT 40.100	MT 100.10 MT 100.50	- disposable - transparent - cyro	- time saving without cleaning - no contamination - facilitated Observation
Multi-use tube	● MMT 40.1	MMT 100.1	- Inside of metal - changeable coupling and sealing	- multi-use
Steril tube	MT 40.10 steril	MT 100.10 steril	- sterilized with gamma radiation	- aseptic condition, i. the absense of living organisms
Disposable tube with titan cutter	MTT 40.10	MTT 100.10	- beater made of Titanium - gray bottom compared to standard tube	- suitable for elemental or heavy metal impurities monitor
Disposable tube with curved cutter	MT 40P.10 MT 40P.100		- Curved beater	- create vortex - effective flow

Application

Milling Chocolate with A11 and Tube Mill

Trial no	Sample+Weight	Device Used	Coolant	Speed	Process Time
1	Silik Oreo (30.1g)	Tube mill 100 + MT 40.10	Dry ice	10000 RPM	30 sec
2	Silik Oreo (33.5g)	A11 mill + A 11.1 blade	Liquid nitrogen	28000 RPM	20 sec



Tube mil results



A11 results


/// Dispersing

Unit operation

/// Mixing and homogenising

› Mixing:

- The process of introducing power into a system, usually for the purpose of producing greater material uniformity



› Homogenising:

- The process of introducing power into a system in order to achieve a specific level of uniformity

Disperse Systeme (Dispersion)

/// Two or more mixable phases

› Classification of disperse systems according to physical state

Dispersion		Disperse phase		
		Solid	Liquid	Gas
Continuous phase	Solid	Dry mixable Media: Granite,	Kneadable media	Porous solid: building material
	Liquid	Suspension: milk	Emulsion: hand cream	Foam: soap foam
	Gas	Aerosol: dust, smoke	Aerosol: mist, spray	Gas mixture (always homogeneous)


I K A U l t r a - T u r r a x - D e v i c e



designed to work perfectly


Rotor-stator system

/// Structure



Rotor-stator system

/// Principle




IKA dispersing technology works by using the rotor-stator principle. The system consists of a rotor within a stationary stator.

Due to the high circumferential speed, the medium to be processed is drawn axially into the dispersion head and then forced radially through the slots in the rotor-stator arrangement.

The high speed and minimal gap between the rotor and stator produces extremely strong shear forces which results in better dispersion.

IKA Dispersers

/// Nomenclature



IKA

designed to work perfectly

Dispersers: Ultraturrax IKA

T10 Basic
0,5-100 ml



T18 Digital
1-1500 ml



T 25 Digital
1-2000 ml



T 50 Digital
0,25-30 litros



IKA+

Scale-up principle
IKA® dispersers have a high degree of flexibility and scalability. Therefore, ensuring reliable scale-up by offering the possibility to work with the same method from formulation development to production.

1:50

T 10 basic | 0,5 – 100 ml | 5000 mPas

T 65 basic | 2 – 50 l | 5000 mPas

UTTD
20-50 ml





designed to work perfectly

Element Disperser Made in Brazil



Ident. Number	Item Description	Product Group	Purchasing Country
3370100	S 10 N - 10 G	Dispersers	IKA BRASIL
3304000	S 10 N - 5 G Dispersing element	Dispersers	IKA BRASIL
004640	S 18 N - 19 G	Dispersers	IKA BRASIL
593400	S 25 N - 18 G	Dispersers	IKA BRASIL
1713800	S 25 N - 25 F Dispersing element	Dispersers	IKA BRASIL
1713300	S 25 N - 25 G	Dispersers	IKA BRASIL
3003900	S 50 N - G 45 F	Dispersers	IKA BRASIL
3003000	S 50 N - G 45 G	Dispersers	IKA BRASIL
3003300	S 50 N - G 45 M	Dispersers	IKA BRASIL

T 10 Basic

	T 10 basic T 10 standard		
Dispersing element	S 10 N – 5 G	S 10 N – 8 G	S 10 N – 10 G
Ident. No.	0003304000	0003305500	0003370100
Working range	0.5 – 10 ml	1 – 50 ml	1 – 100 ml
Stator diameter	5 mm	8 mm	10 mm
Rotor diameter	3.8 mm	6.1 mm	7.6 mm
Gap between rotor and stator	0.1 mm	0.25 mm	0.2 mm
Min. / max. immersion depth	20 / 75 mm	20 / 95 mm	20 / 100 mm
Shaft length	92 mm	115 mm	115 mm
Materials in contact with medium	PTFE, AISI 316L	PTFE, AISI 316L	PTFE, AISI 316L
pH range	2 – 13	2 – 13	2 – 13
Suitable for solvents	yes	yes	yes
Max. temperature	180 °C	180 °C	180 °C
Sterilization methods	all methods	all methods	all methods

①


②

③



T 18 Digital

T 18 digital		
Dispersing element	S 18 N – 10 G	S 18 N – 19 G
Ident. No.	L004639	L004640
Working range	1 – 100 ml	10 – 1500 ml
Stator diameter	10 mm	19 mm
Rotor diameter	7.5 mm	12.7 mm
Gap between rotor and stator	0.35 mm	0.4 mm
Min. / max. immersion depth	25 / 70 mm	35 / 170 mm
Shaft length	108 mm	204 mm
Materials in contact with medium	PTFE, AISI 316L	PTFE, AISI 316L
pH range	2 – 13	2 – 13
Suitable for solvents	yes	yes
Max. temperature	180 °C	180 °C
Sterilization methods	all methods	all methods




T 25 Digital

	T 25 digital			T 25 digital		
Dispersing element	S 25 N – 8 G	S 25 N – 10 G	S 25 N – 18 G	S 25 N – 25 G	S 25 KV – 25 G	S 25 N – 25 F
Ident. No.	0001024200	0000594000	0000593400	0001713300	0002466900	0001713800
Working range	1 – 50 ml	1 – 100 ml	10 – 1500 ml	50 – 2000 ml	50 – 2000 ml	100 – 2000 ml
Stator diameter	8 mm	10 mm	18 mm	25 mm	25 mm	25 mm
Rotor diameter	6.1 mm	7.5 mm	12.7 mm	17 mm	17 mm	18 mm
Gap between rotor and stator	0.25 mm	0.35 mm	0.3 mm	0.5 mm	0.5 mm	0.5 mm
Min. / max. immersion depth	27 / 85 mm	22 / 85 mm	40 / 165 mm	40 / 165 mm	40 / 225 mm	40 / 165 mm
Shaft length	108 mm	105 mm	194 mm	194 mm	270 mm	194 mm
Materials in contact with medium	PTFE, AISI 316L	PTFE, AISI 316L	PTFE, AISI 316L	PTFE, AISI 316L	FFPM / SiC, AISI 316L	PTFE, AISI 316L
pH range	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13
Suitable for solvents	yes	yes	yes	yes	yes	yes
Max. temperature	180 °C	180 °C	180 °C	180 °C	220 °C	180 °C
Sterilization methods	all methods	all methods	all methods	all methods	wet chemical	all methods



T 50 Digital







T 50 digital

Dispersing element	S 50 N – G 45 G	S 50 N – G 45 M	S 50 N – G 45 F
Ident. No.	0008003000	0008003300	0008003900
Working range	0.5 – 20 l	0.5 – 15 l	0.25 – 10 l
Stator diameter	45 mm	45 mm	45 mm
Rotor diameter	36 mm	40.5 mm	40 mm
Gap between rotor and stator	0.5 mm	0.25 mm	0.5 mm
Min. / max. immersion depth	70 / 250 mm	70 / 250 mm	70 / 250 mm
Shaft length	300 mm	290 mm	290 mm
Materials in contact with medium	PTFE, AISI 316L	PTFE, AISI 316L	PTFE, AISI 316L
pH range	2 – 13	2 – 13	2 – 13
Suitable for solvents	yes	yes	yes
Max. temperature	180 °C	180 °C	180 °C
Sterilization methods	all methods	all methods	all methods



designed to work perfectly

PN	Descrição	Chave plana (25001829) 	Chave do eixo (25001830) 	Chave plana para geradores finos (25005638) 	*NOVO* Silentstream (3754000) 
0593400	S 25 N - 18 G	X	X		X
0594000	S 25 N - 10 G	X	X		
1024200	S 25 N - 8 G	X	X		
1713300	S 25 N - 25 G	X	X		X
1713800	S 25 N - 25 F		X	X	X
2348000	S 25 KV - 18 G	X			X
2404000	S 25 KV - 25 F		X	X	X
2466900	S 25 KV - 25 G	X			X
2563000	S 25 KV - 25 G-IL	X			
2830200	S 25 KV - 25 F-IL	X	X	X	
3304000	S 10 N - 5 G	X	X		
3305500	S 10 N - 8 G	X	X		
3370100	S 10 N - 10 G	X	X		
4446500	S 10 N - 8 G - ST		X		
4446700	S 10 N - 10 G - ST	X	X		
4446900	S 25 N - 8 G - ST	X	X		
4447100	S 25 N - 10 G - ST	X	X		
4447300	S 25 N - 18 G - ST	X	X		X
4447500	S 25 N - 25 G - ST	X	X		X
20008825	S 25 KV 2802	X			
20015397	S 25 KV - 26 LR	X			

/// Separation:
Rotary Evaporator made in Brazil

IKA

designed to work perfectly

IKA Rotary Evaporators Overview



Used for manual and standard distilling applications.

/// **The Beginner: RV 3 V** with the new 4 liter heating bath for temperatures of up to 100 °C



Combining both simplicity and steady support through the lift mechanism.

/// **The Athlete: RV 8** with the new 4 l heating bath for water and oil application for temperatures of up to 180 °C




The high-end full automatically controlled version with the 3 liter heating bath for temperatures of up to 180 °C.

/// **The Pioneer: RV 10 auto pro V** with HB 10 heating bath and IKA Vacstar digital pump

Fundamentos de Destilación

El propósito de una destilación es separar materiales diferentes dentro de una cierta mezcla basándose en sus respectivas volatilidades, (puntos de ebullición), a través del proceso de evaporación y condensación (**líquido** para **gaseoso** y **gaseoso** de vuelta al **líquido**).



Fundamentos de Destilación



1. Mezcla a ser destilada

2. Mezcla colectada después de la condensación

3. Condensador dónde ocurre la circulación del agua refrigerante del sistema- Condensación

4. Entrada del vacío (presión negativa)

5. Adición automática de producto (opcional)

6. Rotación del frasco de evaporación



designed to work perfectly

Fundamentos de Destilación

Por que necesitamos de un baño de calefacción?

Evaporación es generada en parte por aumento de temperatura.

Que efecto tiene el vacío en el sistema y por que lo necesitamos?

Cuanto mayor el vacío más reducimos la temperatura necesaria del baño de calefacción para generar la destilación, que por su vez, es acelerada. El vacío también reduce las chances de dañar algunos materiales sensibles al calor.

Que efecto la velocidad de la rotación del frasco de destilación tiene sobre el experimento?

Esa función aumenta la superficie de contacto de la muestra, de esa forma, aumentando la velocidad y eficiencia de la evaporación.

RV 3 Eco

El Iniciante

- Baño de calentamiento de 4 Lts hasta 100°C
- Operación ambidiestra y con una sola mano del elevador manual.
- Velocidad a escala con potenciómetro

Características técnicas

- Desplazamiento: 150 mm
- Velocidad: 20-300 rpm
- Temperatura de calentamiento: temp. ambiental - 99 °C
- Superficie de Condensación: 1500 cm³




Puntos Importantes para recordar:

Target Group:

- Investigación y desarrollo-Farma
- Laboratorios independientes
- Producción e Ingeniería
- Control de calidad

Principales Aplicaciones:

- Química
- Farmacéutica
- Biotecnología
- Académico



✗ Sistema de control de vacío integrado?

✗ Protección extra o trampa para solventes?




RV 3.4 Woulff bottle

IKA

designed to work perfectly

Conexiones RV 3



Vacío- bomba





designed to work perfectly

RV 10 auto y RV 10 auto Pro

El Buque Insignia

- Baño de calentamiento de 4 Lts hasta 180°C
- Operación intuitiva con menú de navegación
- Procesos automáticos- elevador electrónico

Características técnicas

- Desplazamiento: 140 mm
- Velocidad: 5-300 rpm
- Superficie de Condensación: 1500 cm³
- Interfaces: RS 232-USB
- Rotación horaria y antihoraria
- Parada de seguridad integrada



IKA

designed to work perfectly


Puntos Importantes para recordar:

+ Sistema de medición de temperatura?



+ Biblioteca de Solventes

+ Control de vacío integrado?

+ Display y rampas?




Solvents	
Acetic acid	Bath Temperature (°C) 60
Acetone	Rotation (rpm) 100
Benzene	Vacuum (mbar) 556
n-Butanol	Efficiency (kg) 80
Carbontetrachlor	C-Value (L/kg/K) 2100
	Evaporation Enthalpy (J/kg) 553
	Density (kg/m³) 784
	Acetone
SET = Save ESC = Cancel	



Versiones Intermediarias

RV 8

- Equipado con una elevación semi-automática fácil- Función elevación-extracción de seguridad integrada.
- Control óptimo a través de las pantallas digitales
- Funcionamiento sencillo para el usuario y destilación del disolvente necesario con exactitud.
- En todos los paquetes RV 8 se incluye una botella Woulff.



Versiones Intermediarias

RV 10 Digital

- Baño de calentamiento de agua / aceite universal
- Elevador motorizado con función de parada de emergencia
- Interfaz RS 232
- Volumen del baño de calentamiento optimizado para un calentamiento rápido
- Arranque suave
- Rotación a Intervalo a izquierda-derecha
- Función de temporizador
- Circuitos de temperatura de seguridad
- Mecanismo de liberación para soltar la cristalería fija





designed to work perfectly

Aplicaciones

1. Reciclaje de Solventes comunes en el laboratorio: puede fácilmente ser empleado para extraer acetona pura de la acetona utilizada en laboratorio. El circuito de seguridad ajustable, las tapas y alzas de seguridad son características esenciales para dar el soporte necesario en el trabajo da reciclaje de disolventes.

2. Extracciones: aceite de planta o cualquier sustancia semejante. No importa cual el **tamaño** de la muestra, el condensador para reflujo asegura los mejores resultados.

Retención de sobras de metanol na fabricação de Biodiesel

3. Concentración de activos: Elevando la concentración de agentes médicos como aceite de eucalipto es muy seguro y fácil cuando se trabaja con el evaporador RV 10. Especialmente cuando se trata de muestras caras y debe evitarse los tan comunes “impactos”. Con el nuevo sistema RV 10 de desconexión automática y mecanismo de elevación del frasco de evaporación aseguran que impactos no ocurran.

4. Secado de polvos (Principalmente en la área farmacéutica): intervalos entre rotación horaria y anti-horaria. Así es optimizado el secado de polvo para todos los tipos de pigmentos y polvos considerando su frasco específico.

5. Separación de mezclas no azeotrópicas (que no se separan): agua y acetona, Algunas de las características únicas del RV 10 son la alta precisión de las configuraciones de temperatura del.

6. Destilación de sustancias sensibles a temperatura a través del vacío: Con la precisión del RV 10 y constante control del sensor de temperatura del baño, la destilación de sustancias sensibles a temperatura es fácilmente realizada a través del vacío. Un ejemplo es la destilación de ajeno, lo cual, por medio del vacío, conserva el aroma.



designed to work perfectly

Comparación

[PC_RV_20181212-EN.xlsx](#)

IKA

designed to work perfectly

IKA: TODO LO QUE NECESITAS EN UN SOLO LUGAR





designed to work perfectly

Gracias!
Obrigada!
Thank you!

Dudas, preguntas?
Escribanos: sales@ika.net.br